



1441 Wolf Creek Trail  
P.O. Box 302  
Sharon Center, Ohio 44274-0302

PHONE: 1-866-860-8447  
FAX: 1-330-239-4465  
www.transmotionmedical.com

## NEWS RELEASE

October 27, 2009

For more information:

Betsy Mikalacki, 866-860-8447, Ext. 177  
betsymikalacki@transmotionmedical.com

### **TransMotion Medical's TM9 – New Ophthalmology Package**

**SHARON CENTER, OH** – TransMotion Medical, a leading U.S. manufacturer of mobile, motorized stretcher-chairs for hospitals and outpatient surgery centers introduces a new package for ophthalmic surgeries. The **TM9 Ophthalmic Procedure Chair Package** combines the benefits of the TMM5-TB surgical stretcher-chair with options and accessories designed specifically for ophthalmic procedures. The exceptional advantage – the package costs less, up to \$300 less, than ordering the chair with these options and accessories individually.

The TM9 Ophthalmic Procedure Chair Package includes:

- **TMM5-TB Stretcher-Chair** (500 lb. patient capacity) features battery powered, motorized positioning for a low 23" egress for patients, no transferring or manual positioning of the patient by staff, and a high 39" height for superior surgeon access
- **Adjustable Slope Wrist Rest** for surgeon comfort
- **Air Delivery** for the surgical procedure
- **Belt** to keep the patient securely in place during transport and procedures
- **Telescoping I.V. Pole** with 4 insert wells on the frame and can clip under the chair for convenient storage
- **O<sub>2</sub> Tank Holder** for both D and E size oxygen tanks
- **Surgical Support Leg** allowing added stability during procedures

"In addition to this convenient package, we believe the TM9's small footprint will make a large impression on the Ophthalmic market," said Ryan Yerkey, C.O.O. of TransMotion Medical. In its upright position, the chair requires less space than a traditional stretcher and it is easy to navigate in close quarters. The concept is One Patient, One Surface™ -- from pre-op, to procedure, to post-op. By reducing the number of transfers, it streamlines the patient process throughput for increased productivity and reduces patient fall and staff injury risk for a return on investment that is measured in months not years." For more information, call 1-866-860-8447 or visit [www.TransMotionMedical.com](http://www.TransMotionMedical.com).

*TransMotion Medical, Inc. manufactures the most advanced selection of motorized medical stretcher-chairs on the market. Available in more than 30 versatile configurations, this mobile line converts from chair to procedure/stretcher at the push of a button. Designed to reduce, if not eliminate, the need to transfer patients. One patient, one surface™ -- from admission to discharge. TransMotion Medical stretcher-chairs can be found in the top-rated hospitals and surgical centers across the United States as well as international markets. Our mission is to develop for our customers the most innovative, high quality and value specialty treatment chair for safe patient handling. Visit our web site [www.TransMotionMedical.com](http://www.TransMotionMedical.com) to learn more about our products, request a free demonstration or to obtain a quote.*

