



SURE-CHEK® FUSHION II FABRIC WITH 4-WAY STRETCH

Available Colors:



Standard/Default Color



Sure-Chek® Fushion II brand medical fabrics are the original healthcare fabrics that are patented for use as protective covering fabrics and today are the standard in the industry. Sure-Chek® medical fabrics were developed in cooperation with healthcare professionals to provide fabrics that are comfortable, safe and effective for use in healthcare facilities. **These thermoplastic fabric composites are strong, durable, antimicrobial, antistatic, flame resistant, stain, odor, and fluid resistant.**

This healthcare fabric complements pressure management cushions to provide a comfort enhanced product. The protective polymer membrane and the highly stretchable fabric work together to offer protection and comfort.

Changes in demands of the healthcare market have required ever softer, stretchable fabrics. These fabrics perform best on products that are designed to achieve low interface pressure readings. The polymer membrane accommodates the high stretch and provides the recovery properties of this fabric. The stretching feature enables patient immersion into the support surface. The result is the inhibition of pressure sore formation. It is this synergy between the fabric and the cushion that gives the desired outcome.

Sure-Chek® brand medical fabrics are enhanced by a controlled release system, which enables the “slow release” of a specially formulated antimicrobial system to the fabric surface, lasting for the life of the medical fabric.

SOILS OR STAINS: Use neutral soap suds and lukewarm water. Rinse with water and allow to dry.

Do not use harsh cleansers, bleach, solvents or detergents. May also Dry Clean or Autoclave.

DISINFECTION: Dilute disinfectants and/or germicides as specified on the product label. Do not use in concentrations that exceed those that are recommended by the manufacturer. See complete Cleaning Instructions for list of recommended disinfectants.

WWW.TRANSMOTIONMEDICAL.COM

866-860-8447 • info@transmotionmedical.com
Sharon Center, OH