




# SEAT TILT STRETCHER-CHAIRS

- TransMotion Medical manufactures the ONLY STRETCHER-CHAIR with the power to RE-ADJUST PRESSURE POINTS under the patient.

- In 2007, the Center for Medicare and Medicaid Services (CMS) reported 257,412 cases of preventable pressure ulcers as secondary diagnosis. CMS estimated the average cost at \$43,180 per hospital stay.

- As of Oct. 1, 2008, CMS began the Hospital Acquired Condition (HAC) initiative. Medicare will NO LONGER PAY FOR DECUBITUS ULCERS ACQUIRED DURING A PATIENT'S HOSPITAL STAY.










**QUESTION:** Is your hospital doing everything necessary to prevent decubitus ulcers?

**SOLUTION:** TransMotion Medical TMM4-X Series Multi-Purpose Stretcher-Chair or the TMM5-X Series Stretcher-Chair with Articulating headpiece – both featuring the Motorized Seat Tilt Option. Perfect for the in-patient who is required to get out of bed per doctors orders. Smooth lateral transfers and motorized positioning to adjust pressure points.

- › INPATIENT FLOORS
- › MEDICAL-SURGICAL UNITS
- › ICU/CCU
- › SKILLED NURSING FACILITIES
- › SPINAL CORD INJURY
- › GERIATRIC

# SUBJECT #1: FEMALE, 5'2", 123 LBS.

POSITION	PHOTO	PRESSURE MAPPING
<p><b>B – SUPINE POSITION</b> Position in which lateral transfer is easily accomplished</p>		
<p><b>Z – STANDARD SEATED POSITION</b> Typical seated position for feeding, physical therapy. Pressure map shows pressure centralized on the sacral bone.</p>		
<p><b>O – SEAT TILT POSITION</b> By utilizing TransMotion Medical's unique POWERED SEAT TILT. Pressure is significantly reduced to the sacral bone area. By adjusting the tilt of the seat by 1", if significantly shifts pressure points, and relieves sacral pressure.</p>		



X PENDANT

